

ANALIZA BIAŁEK



Scan for verify



100018025

## Laboratory of Analysis

**Client:** Anonymous  
**Order number:** 100018025  
**Name:** SELLIZA TESTO PLUS  
**Number of samples:** 1 ampoule

**Batch number:** 0S25P1543  
**Receive date:** 2026-09-03  
**Start date:** 2026-09-11  
**End date:** 2026-09-13

## Sample identification

| Substance | Method of taking a sample   | Form                                  |
|-----------|-----------------------------|---------------------------------------|
| -----     | Sample provided by customer | Transparent liquid in a glass ampoule |

## Results

| Substance                     | Sample name              | Content     | Purity |
|-------------------------------|--------------------------|-------------|--------|
| Testosterone Propionate       | SELLIZA TESTO PLUS 400MG | 73,18mg/ml  | -----  |
| Testosterone Phenylpropionate | SELLIZA TESTO PLUS 400MG | 74,57mg/ml  | -----  |
| Testosterone Cypionate        | SELLIZA TESTO PLUS 400MG | 96,95mg/ml  | -----  |
| Testosterone Enanthate        | SELLIZA TESTO PLUS 400MG | 146,84mg/ml | -----  |

## Research comment

The test sample was applied in the amount of 3µl to the BioBasic 8 Thermo Scientific column. During the course of the analysis with the method a linear gradient in the solutions of H<sub>2</sub>O and Acetonitrile with the addition of 0.1% trifluoroacetic acid, peaks of the sought substances appeared. Their contents, when compared to the reference chromatograms, are shown in the table above. Height the measurement error is +-1mg. The actual state of the sample sent before the test is presented on the photo below.

## Photo of the test sample




## Password

IY73IJY3

**Analiza Bialek sp. z o.o.**  
**KRS:** 0000966168  
**NIP:** 8992922970  
**REGON:** 521751875

**Headquarters and Laboratory:**  
ul. Duńska 9, 54-427 Wrocław  
Budynek DELTA

**Lab technician:**  
Analysis and research specialist  
mgr Arkadiusz Zajac



ANALIZA BIAŁEK




Scan for verify



100018025

## Chromatogram



## Method

| Time  | H2O + 0,1% TFA | ACN + 0,1% TFA |
|-------|----------------|----------------|
| 0min  | 60%            | 40%            |
| 2min  | 60%            | 40%            |
| 27min | 5%             | 95%            |
| 30min | 5%             | 95%            |

## Password

IY73IJY3

**Analiza Bialek sp. z o.o.**

**KRS:** 0000966168

**NIP:** 8992922970

**REGON:** 521751875

**Headquarters and Laboratory:**

ul. Duńska 9, 54-427 Wrocław

Budynek DELTA

**Lab technician:**

Analysis and research specialist

mgr Arkadiusz Zajęc